**That's comfort!**

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**Apply and remove your compression stockings in 5 easy steps!**

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**Instructions for use**

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**Mode d'emploi**

**The one and only tool for putting on and taking off all types of compression stockings**

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# APPLY / ANZIEHEN / ENFILER

1

Place Steve in front of you on a coffee table and stand behind Steve. Position the stocking above a Steve with the heel towards you.

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2

Fold the top edge of the stocking over the frame. Strip the stocking up until at least the heel reaches the top. Note: the exact point differs per user.

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3

Sit down on a chair. Place the handle against your knee or your groin (if you cannot stretch your foot). Then put your whole foot into the stocking.

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4

Roll Steve towards you so that it is standing flat on the ground under your knee. Then push your heel on the ground(1.). Then pull Steve slowly upwards with up and down movements (2.). Don't overstretch the sock.

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5

Put Steve to one side and rub the stocking into the right place.

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1

Slide the ring with the opening facing away from you with the readable text between the stockings and your leg.

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2

Slide the ring down to the bottom of your calf. Make sure the ring is not hanging crooked.

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3

Place Steve straight onto the ring and push it down.

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4

You can roll Steve away from you as soon as the wheels touch ground.

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5

Keep pushing Steve away from you and pull your leg out of the stocking.

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